



# iSyncBrain<sup>®</sup>

# Advanced EEG Analysis Platform

A.I. automated QEEG analysis platform employing age, sex-classified normative Database.



Cloud Service Access anywhere, anytime



Fast & Easy Upload, receive report in minutes



Science-based Advanced analysis & provision of references



### 🚯 Specialty of our QEEG analysis platform, iSyncBrain



#### DisyncBrain Protocol

- National standard reference data project has been conducted for 8 years at Seoul National University's EEG Data Center.
- ✓ The analysis results report of the EEG file can be exported as a PDF file.

10 Min				
EEG test				
EEG is measured for over 3 minutes in the awake resting state. (Eyes closed / Eyes open)				







iSyncBrain Series



# (1) Age, sex specified EEG standard reference DB production pipeline



# Dirmative data base

- ✓ Standard reference database project has been conducted for 8 years at Seoul National University's EEG Data Center.
- ✓ We have built the world's largest QEEG healthy normative database.
- ✓ We provide EEG analysis results in comparison with the age and sex (4~82 years old) classified QEEG Norm DB.

	Male	Female	Total
Child & Adolescent(4-18)	302	315	617
Adult(19-84)	251	421	672
Total	553	736	1,289

### A.I Automatic denoising

- A.I automatically removes noise components caused by the eye movement, muscle tension, electrocardiogram and pulse, movement, etc.
- It also dramatically reduces your analysis time through machine learning based noisy epoch rejection and AMICA noise components rejections, which normally is done manually.
- iSyncBrain's automatic noise reduction functionality is patented.
  [Patent No. 10-2077605 'Method and apparatus for an automatic artifact removal of EEG based on a deep learning algorithm' 2020.02.10]
- ✓ You may also manually denoise your raw data through iSyncBrain platform.



M



#### Demponent Data

- Time series, Power spectrum (PSD), Topomap, MNI Source location for each components extracted by AMICA are displayed.
- <Artifact Component No.> displays the components determined as noise components by the AI automatic analysis engine



# Band power – Topomap

- EEG power is quantified by frequency 1 Hz bin (1~45Hz) and band specific power.
- Absolute/relative power z score topomaps are displayed.
- Z scores are color coded after comparisons with age/sex classified norm.



# Dewer spectrum

- Power spectrum results show the frequency spectra of the EEG for each channel in units of power spectral density.
- ✓ For visualization, it can be viewed in linear scale (uV<sup>2</sup>/Hz) or log scale (dB/Hz).
- The difference value between left and right symmetric channels is shown as a pair comparison graph.
  If you click the compare button, you can check if there is an asymmetry between the mirrored channels as shown below.





#### 🚯 Power ratio\_TBR / TBR2 / TAR / DAR Theta/beta2 ratio (TBR2) Theta/beta ratio (TBR) TBR2 (Theta Beta2 Ratio) TBR (Theta Beta Ratio) Resting state Resting state Theta-beta2 ratio Theta-beta ratio Associated with ADHD Associated with memory arousal level encoding in cognitive disorder Theta/alpha ratio (TAR) Delta/alpha ratio (DAR) TAR (Theta Alpha Ratio) DAR (Delta Alpha Ratio) Resting state Resting state theta alpha ratio. Delta to Alpha ratio Related to cognition, Related with cognitive deficiency and daily behavioral learning and memory in MCI, Alzheimer's, etc. disabilities after stroke 🚯 Source ROI power (sLORETA) & connectivity(iCoh) Absolute Relative Delta Theta Alpha1 Alpha2 Beta1 Beta2 Beta3 Gamma hypo connection hyper connection Each square represents z score, Alpha1 color coded through comparisons of each ROI's sLORETA activity with the norm. erior Cingulate Anterior Cing Cingul Temporal Pc al Pole-rh ior Temporal-rh vior Tempor Middle Temporal-rh Inferior Temporal-rh Middle Tempora The line shows z score color Inferior Tempor coded imaginary coherence Transverse Tempora nsverse Temporal-rh Banks of Superior Temporal Sulcus anks of Superior Temporal Sulcus-rh between the ROIs. iform-rh orhinal-rh hippocampal-rh al-th z-score All results are displayed in each 0.0 frequency bands. -1.31 -0.65 0,0 0,65 1.31

# (1) 3D View Voxel power / Connectivity

- ✓ Voxel power can be visualized in 3D interactive viewer, z score color coded for each frequency bands.
- ✓ Z score color coded connectivity between source ROIs are also visualized in 3D interactive viewer.
- 3D interactive viewer shows source level activation in any outer directions and inner directions, which enables inspection of medial side activations.

[3D Voxel Power(Abs. / Rel.)]



[3D Voxel Power Left/Right side]



[3D Viewer\_Settings]

Нуро

Hyper

Connection

[3D Connectivity]





# (1) EEG Group analysis

- Group statistics provides analysis of differences before and after the treatment for an individual patient within the group, and between the groups,
- Group statistics also provide independent differences between 2 groups, or 3 groups(ANOVA) which may  $\checkmark$ include control group recruited from iSyncBrain normative library.
- Analysis data can be exported and used as statistical data for research papers.  $\checkmark$



### (1) EEG Group analysis Settings













Comparison between two groups using each group's pre-post differences.



iMediSync, Inc. H. www.imedisync.com T. +82-2-747-7422 E. isyncbrain@imedisync.com

\*This product is for research and analysis purposes only and is not used for diagnostic purposes.

\*This print is a literary work protected by the Copyright law. Reproduction, public transmission, distribution, translation, or conversion to electronic media or machine-readable form of part or all without the prior written consent of iMediSync Co., Ltd. is prohibited.