


# QEEG-based discernment of dementia pathologies through machine learning: Lewy body and Alzheimer's disease

P2-140



Taegyun Jeong<sup>1</sup> (MSc), Ukeob Park<sup>1</sup> (MSc), Seung Wan Kang<sup>1,2</sup> (MD, PhD)

(1) iMediSync, Inc., Seoul, South Korea, (2) Data Center for Korean EEG, College of Nursing, Seoul National University, Seoul, South Korea

*“Identify dementia pathologies through QEEG”*

## INTRODUCTION

- **Alzheimer's disease (AD)**, the most common cause of dementia, destroys nerve cells in the brain through accumulation of **beta amyloid plaques**.
- However, various other pathologies of dementia also exist, such as accumulation of **Lewy body (LB) peptides**.
- Different dementia pathologies carry distinguishing symptoms and treatment methods; hence we cannot disregard the importance of correct identification of the pathology.
- Positron emission tomography (PET) is widely used to screen for pathologies. However, it is **expensive** and results in **exposure to harmful ionizing radiation**.
- Therefore, the present study proposes a quantitative electroencephalography (**QEEG**)-based screening method which overcomes the disadvantages of current screening methods.




Figure 1. Preprocessing/ denoising procedure

- sLORETA was used to estimate source cortical activity at 68 regions defined by the Desikan-Killiany atlas.
- The final dataset: **N = 104; 30 AD; 74 LB.**
- 20% of the data were randomly selected as test data.
- Tree-based algorithms – Random forest, XGBoost, LightGBM were trained with various sets of key hyperparameters.
- **Feature reduction criteria:**
  - Features with p-value < 0.05.
  - Feature importance determined through Shapley values.
  - Exclusion of gamma band (30-45Hz) features due to its vulnerability to high-frequency noise.

## RESULTS

- The best classification performance was achieved by an XGBoost model, with test accuracy of **85.7%**, **AD dementia (ADD)** sensitivity of **83.3%** and **LB Dementia (LBD)** sensitivity of **86.7%**. 5-fold cross validation accuracy was at **82.1%**.




Figure 2. Top 20 features that mainly affected classification results

Sensor & Source level features		
	True ADD	True LBD
Pred ADD	5	2
Pred LBD	1	13

No. features: 25

Table 1. Confusion matrix

- Majority of the features with high importance values were related to delta (1-4Hz), alpha1 (8-10 Hz), alpha2 (10-12Hz), beta2 (15-20Hz) and beta3 (20-30Hz) frequency bands.
- Through a group comparison, we further verified that **LBD delta and alpha band** (8-12 Hz) powers are **stronger** than that of ADD, and **LBD beta band** (15-30 Hz) power is **weaker** than that of ADD.

## CONCLUSIONS

- The classification model developed in the present study showed a promising classification performance, through quantitative EEG data which is **cheap and harmless to record**.
- Such QEEG-based classification models carry a great potential to replace the PET in future, resolving conventional disadvantages.

## CONTACT

Presenting author: Taegyun Jeong

Email: [jtaegyun@imedisync.com](mailto:jtaegyun@imedisync.com)

Second author: Ukeob Park

Email: [hoanouo@imedisync.com](mailto:hoanouo@imedisync.com)

Corresponding author: Seung Wan Kang

Email: [seungwkang@imedisync.com](mailto:seungwkang@imedisync.com)

**iMediSync, Inc.**

**Website:** <https://www.imedisync.com/en/>

Visit us at booth # **920** if you wish to explore our products or become our research partner.



## METHODS

- EEG data employed in the present study were recorded at electrode locations defined by the international 10-20 system, in eyes-closed resting-state condition.
- Figure 1 summarizes the data preprocessing/ denoising procedure.