QEEG-based differentiation between Alzheimer's disease with or without alpha synucleinopathy

(1) iMediSync, Inc., Seoul, Korea, Republic of (South), (2) Yonsei University College of Medicine, Seoul, Korea, Republic of (South)

INTRODUCTION

• As existence of alpha synucleniopathy could affect progression of cognitive impairment in Alzheimer's disease (AD), it's very important to differentiate comorbidity of alpha synucleniopathy in AD to predict progression.



- Lewy body dementia can be clinically diagnosed through visual hallucination, cognitive fluctuation, and parkinsonism.
- However, there is no biomarker that directly confirms pathological findings, so the sensitivity of diagnosis is very low.
- When Alzheimer's disease and Lewy body dementia coexist, the exacerbation of the disease is accelerated compared to only Alzheimer's disease is occurred.

METHODS

- We gathered 3 type of PET scans [18F-Florbetaben brain amyloid-beta PET, FDG-PET, DAT-PET] for patients complaining of cognitive problems.
- Based on pattern of cognitive impairment and PET scan findings, dementia due to Alzheimer's disease [called pureADD] or Lewy Body (called **pureLBD**) or **mixed** type (AD with alpha synucleniopathy, AD with LB or it is also involved vice versa) were clinically classified.
- We measured 19ch EEG based on international 10-20 system which is on eye-closed resting state.
- To find statistically significance in QEEG, 3 labeled groups were performed for each independent t-test comparison which is provided at iSyncBrain[®]. In each comparison, test of normality for each QEEG features was confirmed.

	PureADD	Mixec		
Age (Mean±std)	75.48 ±8.3	77.13 ±8.32		
Male/Female	10/22	63/90		
Table 1. Participants	information			

PureLBD

78.12

±6.27

41/39

Ukeob Park¹, Seung Wan Kang¹ (MD, PhD), Kyoungwon Baik², MD, Byoung Seok Ye², (MD, PhD), "Identify dementia pathologies through QEEG"



Figure 1. Group analysis between PureADD(G1) and PureLBD(G2); (A) Independent t-test with Power Spectral Density(PSD), (B) Occipital Alpha peak frequency comparison



- alpha with low amplitude.
- theta enhancement.

Figure 2. Topomap Comparison between PureADD(G1) and PureLBD

Unit : μV2	Theta			Alpha1		
	Pure ADD	Pure LBD	P-value	Pure ADD	Pure LBD	P-value
Frontal(mean±std)	5.08±5.64	12.06±11.46	0.0000	2.83±2.03	5.69±6.38	0.0005
Central(mean±std)	4.54±4.74	11.22±11.61	0.0000	2.96±2.25	5.46±5.63	0.0011
Temporal(mean±std)	7.21±8.40	16.55±16.89	0.0002	5.79±4.90	9.39±12.29	0.1120
Parietal(mean± s t d)	5.17±6.72	13.63±15.52	0.0001	3.49±3.10	6.96±8.49	0.0267
Occipital(mean±std)	7.90±8.87	22.12±23.23	0.0000	8.823±8.75	20.61±29.93	0.0018

Table 2. Independent t-test table of frequency band power(Theta(4~8Hz), Alpha1(8~10Hz)) * All statistical test is satisfies normality

P1-153 **iMediSync**

• Pure ADD showed the characteristics that general EEG slowing with low total power(Figure 1. (A), Figure 2), relative delta enhancement and desynchronized

Pure LBD showed the pattern of slow alpha peak frequency(Figure 1. (B) with intact synchronization, relatively higher EEG total power and frontotemporal

Compared to pure ADD, ADD with LB showed mixed pattern of ADD and LDB.

Conclusion

- These findings imply that ADD could have different EEG oscillation characteristics according to the Lewybody disease.
- Next step, we will develop the machine-learning algorithm to discriminate the Lewybody disease based on specific EEG features and validate it.

CONTACT

Presenting author: Ukeob Park Email: hoanouo@imedisync.com Coauthor: Kyoungwon Baik Email: baikkw@yuhs.ac Coauthor: Byoung Seok Ye Email: romel79@yuhs.ac Corresponding author: Seung Wan Kang Email: seungwkang@imedisync.com iMediSync, Inc.

Website: https://www.imedisync.com/en/

Visit us at booth **# 920** if you wish to explore our products or become our research partner.

ALZHEIMER'S ASSOCIATION ALZHEIMER'S

22 ASSOCIATION INTERNATIONAL CONFERENCE®