

QEEG-based TabNet classifier for dementia pathologies: Alzheimer's disease and Lewy body disease



Hyerin Jeong¹, Ukeob Park¹, Byong Seok Ye, PhD, MD², Seung Wan Kang, PhD, MD^{1,3}

(1) iMediSync, Inc., Seoul, South Korea

(2) Yonsei University College of Medicine, Seoul, South Korea

(3) Data Center for Korean EEG, College of Nursing, Seoul National University, Seoul, South Korea

TabNet classifier of Alzheimer's Disease and Lewy Body Disease using QEEG

INTRODUCTION

- We used a QEEG-based algorithm to differentiate Alzheimer's disease(AD) and Lewybody dementia(LBD).
- Despite current diagnostic measures, a significant number of LBD cases are undetected or misdiagnosed as AD [2].
- [1] attempted a study on the classification of AD and LBD using RandomForest. Further, the objective of [3] was to contrast the EEG characteristics among LBD/AD+, 'pure'LBD(LBD/AD-) and AD.
- We aimed to address the diagnostic bias in which approximately 70% of dementia diagnoses are attributed to AD, potentially overlooking the presence of LBD.
- The model's efficacy is further evaluated on mixed dementia cases presenting both Alzheimer's and Lewy body symptoms.

METHODS

- 19 channel EEG data, recorded according to the 10-20 system in a resting state, were processed through iSyncBrain®, and automated QEEG platform, to extract sensor-level features for TabNet structure training.
- To address data imbalance, time window-based augmentation applied to EEG sensor level features, creating a final dataset of 7292 samples (AD = 3376, LBD = 3916).
- The dataset was the divided in an 8:2 ratio for model training and testing,
- Mixed type dementia cases were also used to further validate the model's performance.




Fig 1. Data augmentation on overlapping pureAD and pureLBD

	Train	Valid	Test
Original	96	24	30
Augmentation	5755	1387	-
Total	5851	1411	30

	Train		Valid		Test	
	pureAD	pureLBD	pureAD	pureLBD	pureAD	pureLBD
	2710	3141	656	755	10	20

Table 1. Table of augmented dataset number of Train, Valid and test ; pureAD and pureLBD




Fig 2. The process feature selection is performed at each step using DNN and the optimal hyperparameters




Fig 3. Top10 Global feature importance barplot




Fig 4. Masked Feature subset ; red circle represents the features in fig 3.

RESULTS

- This classifier yielded test accuracy of 80% with pure AD sensitivity at 80% and pure LBD sensitivity 80%.

	Pred pureAD	Pred pureLBD
True pureAD	8	2
True pureLBD	4	16

Table 2. Confusion matrix for the testset of the QEEG TabNet classifier

CONCLUSIONS

- Moreover, the predicted probability for mixed dementia data were prominently concentrated around 0.5.
- This model successfully identifies distinctive markers and differentiates between AD and LBD.

REFERENCES

1. Dauwan, M., van der Zande, J. J., van Dellen, E., Sommer, I. E., Scheltens, P., Lemstra, A. W., & Stam, C. J. (2016). Random forest to differentiate dementia with Lewy bodies from Alzheimer's disease. *Alzheimer's & Dementia: Diagnosis, Assessment & Disease Monitoring*, 4, 99-106.
2. McKeith, I. G., Boeve, B. F., Dickson, D. W., Halliday, G., Taylor, J. P., Weintraub, D., ... & Kosaka, K. (2017). Diagnosis and management of dementia with Lewy bodies: Fourth consensus report of the DLB Consortium. *Neurology*, 89(1), 88-100.
3. Van der Zande, J. J., Gouw, A. A., Van Steenoven, I., Scheltens, P., Stam, C. J., & Lemstra, A. W. (2018). EEG characteristics of dementia with Lewy bodies, Alzheimer's disease and mixed pathology. *Frontiers in Aging Neuroscience*, 10, 190.

CONTACT



jhrin@imedisync.com

seungwankang@imedisync.com

iMediSync Inc.

